

Thyroid Cascader – Quick User Guide

The Thyroid Cascader is an interactive clinical decision-support tool designed to help clinicians navigate thyroid-related diagnostic pathways. It combines a lab-pattern interpreter with structured decision trees for Hypothyroidism, Hyperthyroidism, and Thyroid Nodules, guiding users step-by-step to a diagnostic pattern.

1. Home Screen Overview

The home screen provides two main entry options:

Option 1 – Enter labs

Use this when you already have thyroid function tests (TFTs). Select the status of:

- TSH – Elevated, Normal, or Decreased
- Free T4 – Elevated, Normal, or Decreased
- T3 – Elevated, Normal, Decreased, or Unknown / not measured

Press Interpret to generate a result.

The interpretation will show:

- A diagnostic label (e.g., *Primary hypothyroidism*, *Subclinical hyperthyroidism*, *Euthyroid sick syndrome*, *Central hypothyroidism*)
- A list of corresponding etiologies

If the TFT pattern is consistent with:

- Primary hyperthyroidism (overt), or
- Subclinical hyperthyroidism

an additional button appears:

→ Start hyperthyroid pathway

Pressing this takes you directly into the Hyperthyroid cascade (Option 2).

If the results are physiologically normal or non-specific, an explanatory message appears instead of a diagnosis.

Option 2 – Select the presentation

Use this when beginning with the clinical scenario. Three presentation tiles are available:

- Hypothyroid symptoms
- Hyperthyroid symptoms
- Thyroid nodule

Selecting a tile launches the corresponding stepwise decision tree.

Each path contains one decision per screen and guides you toward likely diagnostic patterns.

2. Navigating a Cascade

Each cascade (hypo, hyper, nodule) unfolds as a series of decision nodes.

At each step:

1. Read the question or instruction.
2. Choose the response that matches your clinical findings.
3. Automatically advance to the next node.

Navigation controls:

- Back – return to the previous step
 - Home – return to the main menu
 - Restart – restart the cascade from the beginning (available on outcome screens)
-

3. Hyperthyroid Cascade – Expanded Version

The hyperthyroid pathway includes:

Branching by initial labs:

- Overt hyperthyroidism (TSH low, FT4/T3 high)
- Subclinical hyperthyroidism (TSH low, FT4/T3 normal)
- TSH-mediated hyperthyroidism (inappropriately normal/high TSH)
- Hyperthyroidism unlikely (normal/high TSH with normal hormones)

Overt hyperthyroidism details:

- Evaluation of Graves disease (clinical features → TRAb → uptake scan)
- Radioactive iodine uptake patterns
 - Diffuse → Graves
 - Nodular → toxic adenoma or multinodular goiter

- **Low uptake hyperthyroidism**
 - Positive antibodies → early Hashimoto (hashitoxicosis)
 - High thyroglobulin → endogenous hyperthyroidism (e.g., thyroiditis)
 - Low thyroglobulin → exogenous thyroid hormone use

These now match the latest flowchart-based decision structure.

4. Outcome Screens

Each cascade ends in a clear diagnostic statement accompanied by a brief clinical note.

Outcome screens provide two actions:

- Restart this cascade
- Choose a different path (return home)

The lab interpreter (Option 1) behaves similarly. When a hyperthyroid pattern is detected, it provides an additional shortcut to launch the detailed hyperthyroid cascade.

5. Clinical Use Notes

- This tool supports—but does not replace—clinical judgment.
 - Always integrate:
 - Full clinical history and exam
 - Exact numeric lab values and reference ranges
 - Institutional imaging and biopsy protocols
 - Consider the effects of:
 - Biotin supplementation
 - Assay interference
 - Non-thyroidal illness (euthyroid sick syndrome)
-

6. Intended Users

- Residents and fellows
- Hospitalists and general internists
- Endocrinology specialists

- **APPs (NP/PAs)**
 - **Clinicians evaluating thyroid-related presentations**
-

7. Disclaimer

This tool is educational and decision-support only. It does not replace clinical judgment, institutional guidelines, or specialist consultation.